

Shoes & Socks

Tips for choosing and using them



Shoes

What to look for



Shopping

- You may want to shop without your child. Take the brace with you and try fitting it into the shoe. You can often buy shoes, take them home for a relaxed fitting session, and return or exchange them if needed.
- For online shopping, consider ordering several brands/sizes, then return the shoes that don't work out.
- Once you find a brand you like, go to that specific brand again for a dependable fit.

Trying the shoe on

- Remove the insole layer that comes in the shoe.
- Find the smallest shoe that can hold the brace. You may have to wiggle and push the brace into the shoe before the heel drops in. Use a fairly aggressive push to get the brace down into the toe box. This extra work means the shoe will be only slightly longer than normal. If the brace slides into the shoe too easily, the shoe may be too large and your child may trip over the toes.
- A shoe horn helps get the braced foot into the shoe.

Brands to consider:

• Retail shoes

Stores offering a wide shoe selection often carry styles that are functionally appropriate, stylish, and less expensive than orthopedic shoes.

- | | |
|-----------|-------------|
| Adidas | New Balance |
| BILLY | Nike |
| Champion | Pediped |
| Converse | PLAE |
| DCs | See Kai Run |
| Dinosoles | Skechers |
| Keen | Stride Rite |
| Merrell | Tsukihoshi |
| | Vans |

Stores

- | | |
|------------|--|
| DSW | ShopKo |
| EMS | Shoes.com |
| Fred Meyer | Target |
| JC Penney | Walmart |
| Kohl's | zappos.com |
| REI | <i>(see their Adaptive page for orthotic-friendly shoes)</i> |

• Orthopedic specialty shoes

These shoe brands are designed specifically to accommodate orthoses:

- Answer2
- Hatchbacks
- Markell Shoe Company



Tip for boots and high-tops

We recommend styles with a side or rear zipper, or some other wide-opening feature.

Questions?

Please contact Customer Support:
800.848.7332 | customersupport@dafo.com

Shoes continued

To help the fit

Consider slightly altering the shoe, especially in canvas shoes with a sewn, overlapped toe box.

- Try snipping a few threads that hold the toe box closed around the attached end of the tongue.
- Look at the shoe as you insert the brace and foot to see where the material stretches as it holds the brace. A minor cut at this point can provide more room.



Only one brace?

If your child wears a brace on only one foot, you and your practitioner can discuss the following options:

- The brace will slightly increase the length (height) of the leg it's on. To maintain even leg height:
 - If the shoes come with a removable flat insole, remove it from the braced side, flip it over, and add it to the un-braced side.
 - Put an orthotic into the un-braced foot's shoe. Our HotDog and PattiBob shoe inserts are ideal for this.
- For a bulkier brace, you may need two different shoe sizes—one for the un-braced foot and a larger size for the braced foot.
- For more room, consider altering the braced foot's toe box as described above.



Socks

What to look for

For optimal comfort, look for socks with these features:

- Fit smoothly against the skin, without bunching or folding
- Wick perspiration away from the skin (cotton/poly blends; CoolMax®)
- Seamless
- Resist wrinkling
- Extra-long for full coverage



Tip for sock compliance

For taller AFO brace styles, consider using a colorful soccer sock (or any tall, seamless, moisture-wicking sock). When worn inside-out, the extra length can be folded down over the brace.



Specialty socks:

These sock brands offer all the recommended features, and are specifically designed to wear with SMO, AFO, and KAFO braces.

- Knit-Rite® SmartKnit® Socks
- AFO Socks
- PediatRx Socks

Non-specialty socks:

Socks that have some or most of the recommended features can offer a great value for the price. A moisture-wicking blend of cotton/poly is key.

- Soccer socks
- Knee-highs
- juDanzky and Target's Cat & Jack have received positive reviews.
- A search on Amazon.com for SMO, AFO, or KAFO socks will bring up a number of options.